

# The Latest and Greatest Restorative Kits

by Dr. Lorin Berland



## CONTEMPORARY CUTTING KIT

Bur Block Kit (0956) - Designed by Dr. Lorin Berland.  
Each kit includes eight Robot Diamond Points.



The *Contemporary Cutting Kit™* is for preparing porcelain and composite restorations.

Dr. Lorin Berland's Contemporary Cutting Kit contains Shofu's Robot Diamond Points and NEW technique cutting diamonds. Use these kits particularly when preparing teeth for porcelain/composite veneers, inlays, onlays, and full jacket crowns. It also functions well for removing old crown and bridge restorations and amalgam.




## CONTEMPORARY FINISHING KIT

Bur Block Kit (0957) - Designed by Dr. Lorin Berland. Each kit includes six Ceramisté Standard, Ultra, and Ultra II polishers, five T & F Hybrid Points, and one Robot Diamond Point.

The *Contemporary Finishing Kit™* is for finishing and polishing porcelain and composite restorations.

Dr. Lorin Berland's Contemporary Finishing Kit contains Shofu's T&F Hybrid Points, NEW technique finishing diamonds and a series of Ceramisté' cups and points in Standard, Ultra and Ultra II grits. Finishing touches are made fast and simple with these uniquely manufactured, durable rotary instruments.

The robotic manufacturing process of all Shofu diamond instruments distributes the precise amounts of natural diamond particles onto a one-piece, stainless steel shank. The optimal proportion between the projected and bonded part of the diamond particle ensures a smooth cutting and long lasting performance, everytime!



The Latest and Greatest Video by Dr. Lorin Berland shows detailed usage for each rotary instrument. In the video, Dr. Berland shows step by step, a full mouth reconstruction from start to finish in just two appointments.

To order your copy of Dr. Berland's restorative kits contact your dealer or SHOFU at 1-800-82-SHOFU. To order your copy of the Latest and Greatest Video, contact Dr. Berland at [lorinberland@mindspring.com](mailto:lorinberland@mindspring.com).